

IMPORTANT INFORMATION FOR


WSDA's Digital Image Pest Identification System DIPIS

Can be found at the below website




<http://agr.wa.gov/PlantsInsects/InsectPests/DigitalImagePestID/default.aspx>

Nikon Scope Digital Sight DS-L2 Website Instructions




Appendix A.	Microscope layout
Appendix B.	DS-L2 Main Unit Layout
Appendix C.	Network and Scope Support

E-mail or instruct viewer to go to the following web address
Yakima: <http://68.116.30.179/> Wenatchee: <http://173.209.172.8/>
Pasco: <http://71.93.21.227/> Olympia: <http://173.10.110.34/>
Brewster: <http://64.139.115.171/>




Click onto“ROOM1”



To enter the “**ROOM**” on the Website the viewer will be required to enter the User Name and Password


User Name: **admin**

Password: **wdda**




Note: **all letters need to be lower case**

Viewer enter the following
User Name: **admin**
password: **wsda**



Note: **all letters need to be lower case**


To view image under scope Click onto ROOM1



ROOM1

ROOM1

To view highest quality image click the circle to left the of "1280*960"



ROOM1

To specify refresh rate of viewed image click the “Interval Timer” Drop down menu





Image on screen
refreshes every
30 seconds


ROOM1

To view magnified image click the box to the left of “E. Zoom”




ROOM1

To move around magnified image click the **blue** box or grey squares



To view saved images Click onto ROOM2


ROOM2



ROOM2

Click onto USB0


USB0



ROOM2

Click Folder to view image files

Folder



ROOM2

Click file to view image

File

Accessed a record medium - Microsoft Internet Explorer

File Edit View Favorites Tools Help

Address <https://fortress.wa.gov/agr/scopetumwater1/cgi-bin/cgi-bin/dispdir?D:\IM090422> Go Links

Google Search Bookmarks Find Sign In

msn Powered by Live Search Search News Entertainment Video Sports 84°F

Nikon Digital Sight DS-L2

Entrance ROOM1 **Reload** ROOM3 ROOM4

CF C: USB0 D: USB1 E: USB2 F:

D:\IM090422			
Name	Date	Time	Comment
D:			
IM090422			
22092359.JPG	2009/04/22	09:24:00	
22092344.JPG	2009/04/22	09:23:44	
22091047.JPG	2009/04/22	09:10:48	


http://yahoo9.shockwave.com Internet

ROOM2


To view image

Click Open *to view image immediately*

Click Save *to the copy image to local drive, shared drive or portable device*



Image



Appendix A. Microscope layout

Ergonomically Designed Stereomicroscope






Figure 5

Appendix B. DS-L2 Main Unit Layout

Backside of DS-L2 main unit



Appendix C. Network and Scope Support

Please contact for Network support: helpdesk@agr.wa.gov (subject line: NIKON)
Helpdesk: (360) 902-2008 or (509) 249-6990

Please call for Scope support: Jenni Cena
360-481-5735


Appendix D. Trouble Shooting

Browser Failure


If, Internet Explorer fails to spool up web address for the L2 controller try Chrome

Browser Troubleshooting

1. Browser freezes or stays in “loading” mode. Have controller turned off and exit browser. Wait 30 seconds. Turn on controller. Wait for image on controller. Launch web browser, again.




2. Only use refresh buttons or reload buttons



Browser Troubleshooting

3. **DO NOT** use back button



4. Only upto 4 devices can be longed into the web browser for each device at one time.